beurer

BY 80

DE	Babywaage Gebrauchsanweisung	2
EN	Baby scale Instructions for use	10
FR	Pèse-bébé Mode d'emploi.....	17
ES	Báscula para bebés Instrucciones de uso.....	24

IT	Bilancia pesaneonati Istruzioni per l'uso.....	32
TR	Bebek terazisi Kullanım kılavuzu.....	40
RU	Весы для младенцев Инструкция по применению	47
PL	Waga niemowlęca Instrukcja obsługi.....	54

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Read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Dear customer,

Thank you for choosing a product from our range. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage, beauty, air and baby.

With kind regards,
Your Beurer team

Getting to know your device

A baby’s physical development requires a great deal of energy. This is why regular weight control is important to put parents’ minds at ease that their baby is getting sufficient nutrition and is developing healthily. With the specially developed baby scales from beurer, you will always be able to keep an eye on your child’s weight to ensure their healthy development.

The baby scale has the following features:

- Curved weighing surface
- Large, easy-to-read LCD display
- Automatic and manual hold function for recording the weight of wriggling babies
- Tare weighing
- 20 kg weight capacity
- 5 g graduation
- Automatic switch-off
- Adjustable between weight units (kg/lb/oz)

Table of contents

1. Included in delivery	11	7. Usage.....	14
2. Signs and symbols	11	8. Cleaning and maintenance	14
3. Warnings and safety notes	11	9. What if there are problems?	15
4. Intended use	12	10. Disposal.....	15
5. Device description.....	13	11. Technical specifications	15
6. Initial use	13	12. Warranty / Service	16

1. Included in delivery

- 1x BY 80 baby scale
- 2x 1.5V AA batteries
- These instructions for use

2. Signs and symbols

	WARNING Warning notice indicating a risk of injury or damage to health
	IMPORTANT Safety note indicating possible damage to the device/accessory
	Product information Note on important information
	Read the instructions
	Manufacturer
	CE labelling This product satisfies the requirements of the applicable European and national directives.
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE
	Do not dispose of batteries containing hazardous substances with household waste
	Dispose of packaging in an environmentally friendly manner
	Certification symbol for products that are exported to the Russian Federation and members of the CIS

3. Warnings and safety notes

WARNING

Safety notes

- Protect the baby scale from knocks, moisture, dust, chemicals, drastic changes in temperature and nearby sources of heat (ovens, heaters).
- Keep packaging material away from children (risk of suffocation).
- Repairs must only be carried out by Beurer Customer Services or authorised retailers.
- Do not press the button violently or with pointed objects.

IMPORTANT

- Do not place any objects on the scale when not in use.
- The precision of the scale may be impaired by strong electromagnetic fields (e.g. induction hobs, mobile phones).
- The scale's maximum capacity is 20 kg (44 lb, 704 oz). The results for weight measurement are displayed in 5 g increments (0.02 lb, 0.20 oz).
- Remove any transport locks.

WARNING

Notes on handling batteries

- If your skin or eyes come into contact with battery fluid, rinse the affected area with water and seek medical assistance.
- Choking hazard! Small children may swallow and choke on batteries. Store batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect batteries from excessive heat.
- Risk of explosion! Do not throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a relatively long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, open or crush the batteries.

4. Intended use

The baby scale is only intended for measuring the weight of babies and small children. The baby scale is not a medical device. The baby scale is only intended for private use. It is not suitable for commercial or clinical use. The baby scale must only be used for the purpose for which it is designed and in the manner specified in the instructions for use. Any form of improper use can be dangerous. The manufacturer is not liable for damage resulting from improper or incorrect use.

5. Device description


1. HOLD button 	4. Weighing surface
2. Display	5. Unit adjuster
3. ON/OFF/TARE button 	6. Battery compartment with battery compartment cover

6. Initial use

Inserting the batteries

Before using the baby scale, you must first insert two 1.5 V alkaline mignon AA batteries into the battery compartment of the baby scale. Proceed as follows:

1. Open the battery compartment cover on the bottom of the baby scale.
2. Insert two 1.5V alkaline mignon AA batteries into the battery compartment. Make sure that the batteries are inserted the right way round (+/-), as shown in the battery compartment.
3. Close the battery compartment cover until you hear and feel it click into place.


Product information

The batteries must be replaced when “L” appears on the display.

Selecting the unit of weight

When supplied to the customer, the baby scale is set to weigh in kilograms “kg”. You can find the unit adjuster on the rear of the baby scale. You can use this to also change the unit of weight to pounds “lb” or ounces “oz”.

7. Usage

1. Place the baby scale on an even, hard surface; a hard surface is crucial for correct measurement.
2. Press the ON/OFF/TARE button  to switch on the baby scale. As soon as “0.000 kg” appears on the display, the baby scale is ready to take measurements.
3. Now place the baby on the weighing surface. Do this very carefully; do not let the baby fall and take particular care with the baby's head.
4. This baby scale has a “HOLD” function (AUTO HOLD and MECHANICAL HOLD). The AUTO HOLD function establishes the baby's weight on the display as soon as he/she is still enough to achieve an exact measured value. If the baby is moving around vigorously, his/her weight cannot be recorded. Therefore, try to calm the wriggling baby to achieve a rapid measurement, but do not touch the baby, as this will affect the weight measurement. As soon as the baby's weight is established on the display, the word “HOLD” stops flashing and the weight is displayed for a further 20 seconds. The scale then switches off automatically.
5. If the baby is wriggling about too vigorously, you can select the MECHANICAL HOLD function by pressing the HOLD button. Within 3 seconds, the word “HOLD” stops flashing and the weight is displayed for a further 20 seconds.

TARE function

If you prefer not to place the baby directly onto the weighing surface, but instead onto a hand towel, for example, you must place the hand towel on the baby scale BEFORE switching it on. If the hand towel and other weights that you would like to remove again during weighing weigh more than 50 g, you can also use the TARE function. Once you have pressed the TARE button, “TARE” is shown in the display. To reset the scale to 0.000 g, press the TARE button after “HOLD” has stopped flashing.

8. Cleaning and maintenance

The baby scale should be cleaned from time to time.

This should be done using a damp cloth and, if necessary, a small amount of detergent.

IMPORTANT

- Never use abrasive solvents or cleaning products!
- Do not submerge the baby scale in water.
- Do not rinse the baby scale under running water.
- Do not clean the baby scale in the dish washer.

If you do not plan to use the baby scale for a long period of time, we recommend that you store it in the original packaging in a dry environment. During storage, do not place any objects on the baby scale. Store the baby scale away from children and pets. Remove the batteries from the baby scale.

9. What if there are problems?

Display	Cause	Solution
Err	1. No flat, stable surface. 2. The maximum weight capacity of 20 kg has been exceeded.	1. Place the baby scale on a flat surface without carpet. Re-tare the baby scale. 2. The maximum weight capacity of 20 kg has been exceeded.
Incorrect weight is displayed.	No flat, stable surface.	Place the baby scale on a flat surface without carpet. Re-tare the baby scale.
Incorrect weight is displayed.	The baby is wriggling too much.	Try to calm the baby so that it stops wriggling.
Incorrect weight is displayed.	Incorrect zero point of scale.	Switch the baby scale off and on again. Wait until the baby scale displays "0.000 kg" and repeat the measurement.
Lo	The batteries in the scale are empty.	Replace the batteries in the baby scale.

10. Disposal

For environmental reasons, do not dispose of the device in the household waste at the end of its useful life.

Dispose of the device at a suitable local collection or recycling point in your country. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

Battery disposal

Empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead

Cd = Battery contains cadmium

Hg = Battery contains mercury

11. Technical specifications

Model:	BY 80
Dimensions:	56cm x 33cm x 5.4cm (W x D x H)
Weight:	approx 1.7kg
Measurement range:	0.05kg to 20kg
Scale interval d:	The display can be read in increments of 0.005kg
Repetition accuracy:	The measuring tolerance for repeated measurements is +/- 1% (several measurements in direct succession on the same scale with position of scale and person the same as far as possible)
Absolute precision:	In comparison to a calibrated weight, the measured value is +/-0.2kg.

12. Warranty / Service

Beurer GmbH, Söflinger Straße 218, 89077 Ulm, Germany (hereinafter referred to as “Beurer”) provides a warranty for this product, subject to the requirements below and to the extent described as follows.

The warranty conditions below shall not affect the seller’s statutory warranty obligations which ensue from the sales agreement with the buyer.

The warranty shall apply without prejudice to any mandatory statutory provisions on liability.

Beurer guarantees the perfect functionality and completeness of this product.

The worldwide warranty period is 5 years, commencing from the purchase of the new, unused product from the seller.

The warranty only applies to products purchased by the buyer as a consumer and used exclusively for personal purposes in the context of domestic use.
German law shall apply.

During the warranty period, should this product prove to be incomplete or defective in functionality in accordance with the following provisions, Beurer shall carry out a repair or a replacement delivery free of charge, in accordance with these warranty conditions.

If the buyer wishes to make a warranty claim, they should approach their local retailer in the first instance: see the attached “International Service” list of service addresses.

The buyer will then receive further information about the processing of the warranty claim, e.g. where they can send the product and what documentation is required.

A warranty claim shall only be considered if the buyer can provide Beurer, or an authorised Beurer partner, with

- a copy of the invoice/purchase receipt, and
- the original product.

The following are explicitly excluded from this warranty:

- deterioration due to normal use or consumption of the product;
- accessories supplied with this product which are worn out or used up through proper use (e.g. batteries, rechargeable batteries, cuffs, seals, electrodes, light sources, attachments and nebuliser accessories);
- products that are used, cleaned, stored or maintained improperly and/or contrary to the provisions of the instructions for use, as well as products that have been opened, repaired or modified by the buyer or by a service centre not authorised by Beurer;
- damage that arises during transport between manufacturer and customer, or between service centre and customer;
- products purchased as seconds or as used goods;
- consequential damage arising from a fault in this product (however, in this case, claims may exist arising from product liability or other compulsory statutory liability provisions).

Repairs or an exchange in full do not extend the warranty period under any circumstances.

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