MicroLife WS 80

- ① ONSET Button
- ② User Button
- ③ Down Button
- ④ Up Button
- ⑤ Battery Compartment
- ⑥ KG / LB / ST Switch

Display

- ⑦ Low Battery Indication
- ⑧ Zero Setting
- ⑨ Weight Reading
- ⑩ User Indication
- ⑪ Gender
- ⑫ Age
- ⑬ Height
- ⑭ Body Fat Reading
- ⑮ Body Water Reading
- ⑯ Muscle Mass Reading
- ⑰ Bone Mass Reading
- ⑱ Body Mass Index (BMI) Reading
- ⑲ Basal Metabolic Rate (BMR) Reading
- ⑳ Error

Read the instructions carefully before using this device.

Dear Customer,
This diagnostic scale allows you to evaluate your body composition. The knowledge of your body composition allows you to combine a correct nutrition and physical activity as essential conditions to improve the quality of life, promote health and prevent pathologies. Detailed user information about our products and about our services can be found at www.microlife.com.
Stay healthy!

1. Important Facts about Body Composition

It is important to know the body composition of an individual in order to be able to estimate his/her health condition. It is not enough to know one's weight but it is necessary to define its quality, which means dividing the human body into its two main components:

- Lean body mass: bones, muscles, viscera, water, organs, blood
- Fat body mass: adipose cells

Body Fat
The body fat mass is calculated by totaling the essential fat (or primary fat) and the spare fat (or storage fat). The essential fat is necessary for the metabolism in some spare fat, and above 70 % is unhealthy. Research shows that 57 % body mass is optimum, but normally there is no problem if the reading is 10 % higher or lower. The body water % is greater in childhood and increases with age. Women, usually having more adipose tissue than men, have a water ratio of about 55-58 % which men have about 60-62 %. Normally below 40 % would be considered too low, and above 70 % too high.

Body water variations that take place in a few hours are completely normal and are due to water retention, since the body water is the only component subject to such rapid variations.

Muscle Mass
Muscle mass is important in determining a healthy body composition. A person with a higher % of muscle mass finds it easier to move, it needs more energy to do so. Exercise is very important in maintaining a healthy body and the muscle mass % is a useful indicator to control it. The normal muscle mass % is between 28-34 % for men and between 28-39 % for women depending on age and physical activity level.

Bone Mass
Bone mass increases rapidly in childhood and reaches its maximum between 30-40 years. It decreases slightly afterwards with increasing age. Healthy nutrition, regular exercise including some level of aerobic building and a good bone mineralization. However, it is otherwise difficult to significantly improve bone mass.

Factors like age, gender, weight, and height have a small influence on the bone mass. Average adult bone mass is 15 % for male and 12 % for female. There are no recognized guidelines or recommendations relating to bone mass.

Body Mass Index (BMI)
The BMI is a simple index of weight-for-height that is commonly used to classify underweight, overweight and obesity. It is calculated by dividing the body weight in kg by the square of the height in m of a person. The WHO BMI Classification is as follows:

- Normal: 18.50-24.99 kg/m²
- Underweight: below 18.50 kg/m²
- Overweight: over 24.99 kg/m²
- Obese: over 29.99 kg/m²

Basal Metabolic Rate (BMR)
The BMR is the amount of energy required by the body when in a state of complete rest to maintain its basic functions. About 70 % of a human's total energy expenditure is due to the basal life processes within the organs of the body. About 20 % of the BMR comes from physical activity and another 10 % from the digestion of food.

Estimation can be acquired through an equation using age, sex, height, and weight. There are no recognized guidelines or recommendations relating to BMR.

2. Suggestions for Use

- Place the scale on a flat, hard surface. Soft, inaccurate flooring (e.g. rugs, carpets, linoleum) can cause the scale to give uneven readings of your weight.
- Step onto the scale with your feet parallel and weight equally distributed.
- Stand still while the scale measures your weight.
- It is important to wear shoes without laces, or to correctly position on the electrodes. Moreover, for the scale to work properly, avoid your legs touching.
- The scale is safe to use as a simple scale so everyone, including children, can use it.
- The precision of the measurement results can not perform a diagnostic measurement:

- pregnant women
- persons on an oral contraceptive medication
- persons fitted with a pace-maker or other forms of electronic, medical equipment
- those with a temperature
- children under the age of 10
- those in dialysis, with oedemas on their legs, people with thrombosis or suffering from osteoporosis

Always weigh yourself on the same scale and at the same time, preferably dressed and before breakfast. To get the best results from your scale, weigh yourself twice, and if the two weights are different from each other, your weight is between 10-100 years and those aged up to 15 minutes after getting up, so that the water can distribute itself around your body.

The surface of the scale can be slippery if wet. If an error occurs in the analysis of your body composition, the scale will turn off automatically. The precision of the measurement results can not perform a diagnostic measurement.

3. Using the Device for the First Time
This scale operates with two replaceable lithium batteries. Before the batteries are used, they must be inserted in the battery compartment. Insert two new batteries in the battery compartment (➔) on the positive pole upwards.

MicroLife WS 80

- ① ONSET Button
- ② User Button
- ③ Down Button
- ④ Up Button
- ⑤ Battery Compartment
- ⑥ KG / LB / ST Switch

Display

- ⑦ Low Battery Indication
- ⑧ Zero Setting
- ⑨ Weight Reading
- ⑩ User Indication
- ⑪ Gender
- ⑫ Age
- ⑬ Height
- ⑭ Body Fat Reading
- ⑮ Body Water Reading
- ⑯ Muscle Mass Reading
- ⑰ Bone Mass Reading
- ⑱ Body Mass Index (BMI) Reading
- ⑲ Basal Metabolic Rate (BMR) Reading
- ⑳ Error

Read the instructions carefully before using this device.

Dear Customer,
This diagnostic scale allows you to evaluate your body composition. The knowledge of your body composition allows you to combine a correct nutrition and physical activity as essential conditions to improve the quality of life, promote health and prevent pathologies. Detailed user information about our products and about our services can be found at www.microlife.com.
Stay healthy!

1. Important Facts about Body Composition

It is important to know the body composition of an individual in order to be able to estimate his/her health condition. It is not enough to know one's weight but it is necessary to define its quality, which means dividing the human body into its two main components:

- Lean body mass: bones, muscles, viscera, water, organs, blood
- Fat body mass: adipose cells

Body Fat
The body fat mass is calculated by totaling the essential fat (or primary fat) and the spare fat (or storage fat). The essential fat is necessary for the metabolism in some spare fat, and above 70 % is unhealthy. Research shows that 57 % body mass is optimum, but normally there is no problem if the reading is 10 % higher or lower. The body water % is greater in childhood and increases with age. Women, usually having more adipose tissue than men, have a water ratio of about 55-58 % which men have about 60-62 %. Normally below 40 % would be considered too low, and above 70 % too high.

Body water variations that take place in a few hours are completely normal and are due to water retention, since the body water is the only component subject to such rapid variations.

Muscle Mass
Muscle mass is important in determining a healthy body composition. A person with a higher % of muscle mass finds it easier to move, it needs more energy to do so. Exercise is very important in maintaining a healthy body and the muscle mass % is a useful indicator to control it. The normal muscle mass % is between 28-34 % for men and between 28-39 % for women depending on age and physical activity level.

Bone Mass
Bone mass increases rapidly in childhood and reaches its maximum between 30-40 years. It decreases slightly afterwards with increasing age. Healthy nutrition, regular exercise including some level of aerobic building and a good bone mineralization. However, it is otherwise difficult to significantly improve bone mass.

Factors like age, gender, weight, and height have a small influence on the bone mass. Average adult bone mass is 15 % for male and 12 % for female. There are no recognized guidelines or recommendations relating to bone mass.

Body Mass Index (BMI)
The BMI is a simple index of weight-for-height that is commonly used to classify underweight, overweight and obesity. It is calculated by dividing the body weight in kg by the square of the height in m of a person. The WHO BMI Classification is as follows:

- Normal: 18.50-24.99 kg/m²
- Underweight: below 18.50 kg/m²
- Overweight: over 24.99 kg/m²
- Obese: over 29.99 kg/m²

Basal Metabolic Rate (BMR)
The BMR is the amount of energy required by the body when in a state of complete rest to maintain its basic functions. About 70 % of a human's total energy expenditure is due to the basal life processes within the organs of the body. About 20 % of the BMR comes from physical activity and another 10 % from the digestion of food.

Estimation can be acquired through an equation using age, sex, height, and weight. There are no recognized guidelines or recommendations relating to BMR.

2. Suggestions for Use

- Place the scale on a flat, hard surface. Soft, inaccurate flooring (e.g. rugs, carpets, linoleum) can cause the scale to give uneven readings of your weight.
- Step onto the scale with your feet parallel and weight equally distributed.
- Stand still while the scale measures your weight.
- It is important to wear shoes without laces, or to correctly position on the electrodes. Moreover, for the scale to work properly, avoid your legs touching.
- The scale is safe to use as a simple scale so everyone, including children, can use it.
- The precision of the measurement results can not perform a diagnostic measurement:

- pregnant women
- persons on an oral contraceptive medication
- persons fitted with a pace-maker or other forms of electronic, medical equipment
- those with a temperature
- children under the age of 10
- those in dialysis, with oedemas on their legs, people with thrombosis or suffering from osteoporosis

Always weigh yourself on the same scale and at the same time, preferably dressed and before breakfast. To get the best results from your scale, weigh yourself twice, and if the two weights are different from each other, your weight is between 10-100 years and those aged up to 15 minutes after getting up, so that the water can distribute itself around your body.

The surface of the scale can be slippery if wet. If an error occurs in the analysis of your body composition, the scale will turn off automatically. The precision of the measurement results can not perform a diagnostic measurement.

3. Using the Device for the First Time
This scale operates with two replaceable lithium batteries. Before the batteries are used, they must be inserted in the battery compartment. Insert two new batteries in the battery compartment (➔) on the positive pole upwards.

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