



# CN04 COMPRESSOR NEBULIZER



It is a compact medical device designed to efficiently deliver physician prescribed medication to the bronchial lung passages. With proper care and use, it will provide you with many years of reliable treatment.

This product is developed for the successful treatment of asthma, allergies and other respiratory disorders. It creates a stream of air that travels through clear tube to the nebulizer. When air enters the nebulizer, it will convert the prescribed medication into aerosol mist for easy inhalation.

Your Compressor Nebulizer should be used under the supervision of a licensed physician and/or a respiratory therapist after reading the instructions manual.

## SPECIFICATION

VOLTAGE	230 VAC $\pm$ %10, 50 Hz
POWER	100 W
MEDICATION CAPACITY	5 ml
PARTICLE SIZE	0.5 - 10 $\mu$ m
MMAD	4 $\mu$ m
NOISE LEVEL	BELOW 60 dBA
MAXIMUM MEDICATION CAPACITY	13 ml (cc)
AVERAGE NEBULIZATION RATE	min 0,2 ml / min
COMPRESSOR PRESSURE RANGE	30 - 36 PSI (max.)
OPERATION PRESSURE RANGE	8 - 16 PSI
LITER FLOW RANGE	6 - 8 lpm
OPERATING TEMPERATURE RANGE	10°C - 40°C
OPERATING HUMIDITY RANGE	10 - 90 % RH
WEIGHT	< 2 kg

## ACCESSORIES

