

10L Oxygen Concentrator OZ-5-01GW0



- PSA (Press Swing Adsorption) Methodology to Deliver Oxygen without any Consumables but Electrical Power
- Delivers a 93%+/-3% Purity Oxygen at the Flow of 10 LPM (liters per minute)
- Goes with a built-in nebulizer.
- Goes with dual outlets, enabled to work as two oxygen concentrator 0-5 LPM at the same time, or work as one oxygen concentrator 0-10 LPM.
- Oxygen Output pressure: 0.05MPa
- Transportation/Storage Humidity: $\leq 80\%$ RH
- Transportation /Storage Temperature: -20°C-60°C

Oxygen Concentration

-Purity of oxygen reaches 93%±3%

Multiple Ways to use

-It comes with two flow-meters having flow from 0-5 liter and 0-10 liter per minute, it can deliver 5LPM oxygen with two persons at the same time, and can deliver 10LPM with one person.
-It also has atomization function

Safety

-The machine has audiovisual alarm for power failure, low oxygen level, low pressure

Product Details



10L Flow-meter

5L Flow-meter



Large LCD display

Product specifications

Dimensions	390*350*720MM
Weight	29.5KG
Flow Specifications	10LPM (Dual 5 LPM each)
Oxygen Concentration	93%±3%
Power Consumption	750 Watts
Sound Level	60dB(A)