

Wrist Blood
Pressure Monitor

UB-533



Wrist Blood
Pressure Monitor

UB-525



Wrist Blood Pressure Monitor UB-533/UB-525

UB-533

All-in-one Wrist Monitor



The UB-533 is designed for 2 users and focuses on greater ease of use with professional quality. It is equipped with functions like CPG, cuff fit and a movement error indicators to support users in accurate measurement.

The %IHB and WHO blood pressure classification indicators provide users with detailed monitoring.

UB-525

Simple Operation
with Accuracy



The UB-525 is equipped with all functions essential for blood pressure monitoring. Its sleek design comfortably fits the wrist and provides accurate measurement.

About "CPG" Function(UB-533)



CPG (Correct Position Guidance) indicates the correct arm position with both the color and position of the light.



Low angle
Incorrect height

Too LOW

Proper angle
Correct height

CORRECT

High angle
Incorrect height

Too HIGH

Specifications

	UB-533	UB-525
Measurement method	Oscillometric measurement	
Display type	Digital LCD, Pressure/Pulse displayed simultaneously	
Measurement range	Pressure : 0 - 299 mmHg Pulse : 40-180 beats / minute	
Measurement accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$	
Power supply	2 x 1.5V alkaline batteries (LR03 or AAA)	
Number of measurements	Approx. 250 measurements, when AAA alkaline batteries are used, with pressure value of 170 mmHg at room temperature of 23°C.	
Wrist circumference	13.5 - 21.5 cm	
Memory storage	60 x 2 users	60
Memory data average	Maximum 60 data average of systolic, diastolic, and pulse values	
Operating conditions	+10°C to +40°C / 15%RH to 85 %RH / 800 hPa to 1060 hPa	
Transport/Storage conditions	-20°C to +60°C / 10%RH to 95%RH	
Dimensions	Approx. 56 [W] x 88 [H] x 21.5 [D] mm	
Weight	Approx. 105g, excluding the batteries	



Correct Position Guidance



Cuff Fit Error Indicator



Movement Error Indicator



%IHB Indicator to show frequency of Irregular Heartbeat



WHO Classification Indicator



Irregular Heartbeat Indicator



2 user memory (60 x 2)



60-memory



Clock & Calendar



Average reading of last 60 measurements

Standard accessories

- Blood pressure monitor
- Batteries for temporary use x 2
- Instruction manual
- Carrying case(UB-533)

AND
A&D Medical

A&D Company, Limited

3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013 Japan
Telephone:[81](3)5391-6132 Fax:[81](3)5391-6148
1-243, Asahi, Kitamoto-shi, Saitama, 364-8585 Japan
Telephone:[81](0485)93-1111 Fax:[81](0485)93-1119
<http://www.aandd.jp>

A&D INSTRUMENTS LIMITED

Unit 24/26 Blacklands Way Abingdon Business Park,
Abingdon, Oxon OX14 1DY UNITED KINGDOM
Telephone:[44](1235) 550420 Fax:[44](1235) 550485
(Authorized Representative Established in the European Community)